



best@buchi No. 19

Nutrition: Analysis of cheese by FT-NIR Technology

This study was conducted to develop a rapid and accurate NIR analysis method for determining fat, moisture, protein and NaCl in cheese. Commercial cheeses (thousand samples) were used to create calibration models. Conventional analysis of cheese ingredients is both time-consuming and expensive. Near infrared (NIR) spectroscopy is fast, inexpensive and can provide useful information on the ingredients of cheese.



Determination of the composition of cheese

Near infrared (NIR) spectroscopy was used to determine the fat, moisture, protein and salt content of hard, semihard and young cheese. NIR spectroscopy is a rapid and precise analysis method and can be used to replace wet chemistry to simultaneously quantify different parameters of cheese. In comparison to traditional wet chemistry, which requires hours for the analysis of each component, the time requirement of this method is slashed to 16 seconds. A set of thousand samples from diverse types of cheese (Emmentaler, Gruyeres, St. Paulin, Raclette, Appenzell, Tisler etc) were used to calculate a calibration for an FT-NIR spectrometer of type BUCHI NIRFlex N-500.

Author: G. Campolongo

Near infrared spectroscopy

Chemical analysis provides a listing of the components in a particular food or ingredient. Given that most of these are organic molecules which absorb light in the near infrared (NIR) spectral region (750-2500nm), the NIR spectrum of a particular food or feedstuff provides information about its chemical composition. Additional knowledge about its physical makeup may also be present. Analytical laboratories have looked to the ability of spectroscopy techniques to characterize the quality and/or authenticity of foods.

Advantages and limits of NIR spectroscopy

The major advantages offered by NIR spectroscopy include the speed of analysis, the fact that no or only little sample preparation is needed (such as mixing or grinding), and the simultaneous measurement of different parameters. On the other hand, NIR is an indirect method which requires calibration. The cost of developing calibrations for complex products is substantial. However, tools for calibration development have been improved over the past five years, and many of the routine problems can now be kept under control (for example Calibration Wizard for automatic calibrations in BUCHI NIRCal). Nevertheless, to obtain a useful and reliable calibration is not a simple task. It is usually a multi-step procedure in which it is unlikely that all variables will be covered by the initial calibration. BUCHI is therefore offering pre-calibrated applications to be used with NIRFlex. The quality analysis in a typical laboratory procedure is limited to a handful of samples a day when classical methods are used. When an error occurs in the production line, the delay may mean that a significant amount of product will be wasted before the results are known and the process

can be corrected. Fast availability of measurement results is therefore highly valuable.

Materials and methodology

1. Properties of interest

Fat, salt, protein and moisture content are essential for determining hard cheese quality in the cheese industry as well as for cheese research. Fat is a very important parameter for the richness of flavor. Moisture, on the other hand, relates to the profitability of the final product, protein and the salt content are essential for the taste and shelf life of the cheese. Therefore, the relevant contents of these ingredients are important quality and monetary parameters.

2. Type of samples

A set of thousand samples from diverse types of cheese (Emmentaler, Tisler, Appenzeller, Cheddar, Gruyère, etc.) were used to calculate a quantitative calibration for moisture, fat, protein and salt.

3. BUCHI NIRFlex N-500 FT-NIR spectrometer

This instrument provides quantitative measurements of food and animal feed samples using light in the 1000 to 2500 nm range. The NIRFlex N-500 is characterized by its excellent wavelength reproducibility, which is crucial for ensuring calibration portability and stability. The instrument is equipped with an extended range in gas detector. It yields spectra up to 1500 data points for each spectrum.

4. Measurement cell NIRFlex Solids

The measurement cell used is the NIRFlex Solids, which measures samples in diffuse reflectance. A rotating sample desk enhances the repeatability of measurement results. Samples are measured in standard glass petri dishes. The sample may also be measured in monouse plastic petri dishes.



5. Software



The software is divided into two separate packages:

For calculation of the calibration the BUCHI Chemometric software NIRCal is used. NIRCal optimizes

the calibration with wizards, 2D and 3D representations and all appropriate algorithms. For cheese calibrations, the PLS (Partial Least Squares) method was

used. Routine measurements are performed with NIRWare, a modular and easy-to-use software package with SOPs for sample analysis, giving step-by-step instructions.

Data collection & development of the calibration

Sample presentation

Samples of about 100 g were grated using a commercial cheese grater and without applying too much mechanical pressure. The samples were transferred to glass petri dishes and measured by means of the NIRFlex Solids measuring cell in diffuse reflectance. A measurement consisted of 64 scans. Each sample was automatically measured three times. The reference values only had to be allocated once per sample.

Calibrations


To establish the calibrations, several different cheese varieties (Tête de Moine, Raccard, Appenzeller, Parmesan, Chester, Gauda, Emmentaler, Edamer, Asiago) were analyzed on the basis of 10-13 different samples each, which had different production dates. The spectra were recorded consecutively using several NIRFlex N-500 units to ensure maximum robustness of the calibration. A common library was created of all the spectra which served as a basis for calculating the calibration. It was found that one calibration is sufficient for each individual parameter, i.e. there is no need to switch between calibrations with high or low content. The samples were measured at room temperature.

Composition of test material

Fat	6 - 37 %
Moisture	27 - 68 %
Protein	4 - 51 %
NaCl	0.1 - 4.0%

Reference methods:

The reference values required for calibration were determined by the following methods: Moisture: drying oven; salt content: titration with silver-nitrate solution; fat: Gerber or Soxhlet method, Protein: Kjeldahl method. The diagram below shows all spectra used for calibration. To achieve a successful calibration calculation, the division into calibration and validation spectra is crucial, beside the correct data pretreatment



The calibrations were calculated using the NIRCal software. Different data pretreatment methods were used in doing so:

Parameter	Method	Range of calibration wavenumber [cm-1]
Fat	PLS	4598 - 10000
Moisture	PLS	4400 - 9000
Fat in dry matter	PLS	4397 - 10000
NaCl	PLS	4230 - 9250

Diagram 1 compares the predicted values with the reference values of the "Fat in cheese" calibration. The strength of the calibration with a good coverage in the range of 12-37% fat and the rather weaker coverage in the lower range are apparent.




Diagram 1: Fat in hard cheese calibration

Diagrams 2 and 3 show the relevant straight calibration lines for the parameters protein and moisture. Both calibrations convince by their good coverage across the entire range.

Results and discussion

The accuracy of the calibrations can be assessed on the basis of the standard error of prediction (SEP) and the regression coefficient of the straight calibration line. The closer this Rcoefficient is to 1, the more accurate is the prediction of new samples or matching with the values of conventional analysis.




Diagram 2: Protein in cheese calibration




Diagram 3: Moisture in cheese calibration

Properties	Number of samples spectra for calibration	Calibration range (%)	SEP	Standard-deviation of reference method	R-Co efficient
Mositure	6253	27 - 68	0.59	0.31	0.996
Fat	5915	6 - 37	0.57	0.28	0.995
Protein	590	4 - 51	0.73	0.35	0.997
Salt (NaCl)	2078	0.1 - 4.0	0.16	0.06	0.973

Table 1: Calibration performance

In addition to the precision presented here, FT-NIR technology also offers numerous practical benefits in routine analytics in the field of quality inspections. They include:

- Versatility: Simultaneous measurements of different parameters
- Rapidity: Results available in a matter of seconds
- Economy: Low operating costs
- Ecology: No chemical waste
- Safety: No chemical reaction
- Simplicity: No specialist training required for operation
- Calibration portability: Fourier Transform Spectrometers allow calibration data to be transferred from one instrument to another

Quality in your hands

Find your BUCHI Contacts:

BÜCHI Labortechnik AG
 CH – 9230 Flawil 1
 T +41 71 394 63 63
 F +41 71 394 65 65
 buchi@buchi.com
 www.buchi.com

BUCHI Italia s.r.l.
 IT – 20090 Assago (MI)
 T +39 02 824 50 11
 F +39 02 57 51 28 55
 italia@buchi.com
 www.buchi.it

BUCHI Hong Kong Ltd.
 HK – Central
 T +852 2389 2772
 T +852 2389 3366
 F +852 2389 2774
 china@buchi.com
 www.buchi.com.cn

Nihon BUCHI K.K.
 JP – Tokyo 110-0008
 T +81 3 3821 4777
 F +81 3 3821 4555
 nihon@buchi.com
 www.nihon-buchi.jp

BUCHI Korea Inc
 KR – Seoul 153-782
 T +82 2 6718 7500
 F +82 2 6718 7599
 korea@buchi.com
 www.buchi.kr

BÜCHI Labortechnik GmbH
 DE – 45127 Essen
 Freecall 0800 414 0 414
 T +49 201 747 490
 F +49 201 747 492 0
 deutschland@buchi.com
 www.buechigmbh.de

BÜCHI Labortechnik GmbH
 Branch Office Benelux
 NL – 3342 GT
 Hendrik-Ido-Ambacht
 T +31 78 684 94 29
 F +31 78 684 94 30
 benelux@buchi.com
 www.buchi.nl

BUCHI Shanghai
 CN – 200052 Shanghai
 T +86 21 6280 3366
 F +86 21 5230 8821
 china@buchi.com
 www.buchi.com.cn

BUCHI India Private Ltd.
 IN – Mumbai 400 055
 T +91 22 667 75400
 F +91 22 667 18986
 india@buchi.com
 www.buchi.in

BUCHI Corporation
 US – New Castle,
 Delaware 19720
 Toll Free: +1 877 692 8244
 T +1 302 652 3000
 F +1 302 652 8777
 us-sales@buchi.com
 www.mybuchi.com

BUCHI Sarl
 FR – 94656 Rungis Cedex
 T +33 1 56 70 62 50
 F +33 1 46 86 00 31
 france@buchi.com
 www.buchi.fr

BUCHI UK Ltd.
 GB – Oldham OL9 9QL
 T +44 161 633 1000
 F +44 161 633 1007
 uk@buchi.com
 www.buchi.co.uk

BUCHI (Thailand) Ltd.
 TH – Bangkok 10600
 T +66 2 862 08 51
 F +66 2 862 08 54
 bacc@buchi.com
 www.buchi.co.th

PT. BUCHI Indonesia
 ID – Tangerang 15321
 T +62 21 537 62 16
 F +62 21 537 62 17
 indonesia@buchi.com
 www.buchi.co.id

BUCHI Russia/CIS
 United Machinery AG
 RU – 127787 Moscow
 T +7 495 981 05 22
 F +7 495 981 05 20
 russia@buchi.com
 www.buchi.ru

We are represented by more than 100 distribution partners worldwide.
 Find your local representative at: www.buchi.com

