

best @buchi



www.buchi.com

Information Bulletin

Number 62/2011

Sample: B S -50mg_11 _00	
Identity:	ok
Expected substance:	014 -51
Found substance:	ok 014 -51
Distance:	ok 0.001021
Spectrum residual:	ok 0.004844
Hitlist	
Hit	Substance
1	014 -51 Topiramate 50mg
2	01400-01638 Topiramate 50mg

NIRFlex
Solids Transmittance

Quality Control of Tablets: Identification of the Active Ingredient Topiramate Using NIR-Spectroscopy

Author:

Rolf Erismann, Cilag AG, CH-Schaffhausen

Scope

The identity of clinical test samples has to be verified regularly (placebo – verum). This task requires not only a fast but above all a nondestructive method. Without doubt FT-NIR transmission spectroscopy is the preferred technology for the analysis of intact tablets and BÜCHI Labortechnik AG has been identified as a competent provider and partner. The NIR experiments led not only to a validated method that could differentiate between the active ingredient and the placebo, but also yielded the prospect of the ability to discriminate between the individual strengths. An emerging task was the need to transfer the calibration from an earlier BUCHI instrument (NIRTab) to the current NIRFlex N-500 model (Fig. 1).

Original Analysis Method

The active ingredient Topiramate is used in coated tablets in strengths of 25, 50, 100 and 200 mg. Essentially Topiramate does not exhibit any UV absorption, such that identification must be undertaken by measuring an IR spectrum. For this test several Topiramate tablets are crushed and dissolved in acetone, centrifuged and subsequently filtered. Applied to NaCl crystals, the IR spectrum is detected and compared with the reference spectrum. The time required for sampling, analysis, cleaning and documentation is 15 minutes and is carried out in the majority of cases for a sequence of eight batches.

NIR-Analysis Method

156 spectra from 36 tablet batches were measured with NIRTab instrument and the calibration prepared. 73 independent spectra were used for the validation (validation set). In the spectra in Fig. 2 it can be seen that there are only minor differences between the four strengths and that the transmission spectra clearly



Fig. 1: NIRFlex N-500 with transmission module




Fig. 2: NIR spectra of Topiramate tablets




Fig. 3: Zoomed pretreated spectra

overlap, even in the subsequently calibrated wavelength range.

Visual differentiation is not possible without further means. A precise and correct interpretation of a NIR spectrum is only possible by the appropriate application of mathematical spectrum preprocessing. These so called pretreatments have a significant impact on the quality of the calibration (e.g. robustness). The chemometric software tool NIRCal offers numerous mathematical pretreatments of the spectra. Using the latest version of NIRCal 5 with Toolbox the analyst can automatically include countless calibrations in the calculation and the selection can be tightly limited based on a quality factor (Q-value). After the measurement of the spectra and the menu based calculation of calibrations, the user must then evaluate the calibrations offered. NIR and chemometric know-how is obtained during a BUCHI training course.

And very important: it is possible to contact the BUCHI application team at any time and obtain substantial support, as one could say, almost a complete calibration. In the example of the Topiramate calibration, the following spectra

pretreatment is used for the cluster method (Fig.3):

1. Smoothing Savitzky-Golay 9 points
2. Normalization between 0 and 1
3. MSC amplification.




Fig. 4: Cluster calibration

The "cluster plot" (Fig. 4) shows that all four tablet strengths are separated. The Q-value at 0.91 is good. All 73 validation spectra have been correctly allocated (correct prediction). The BUCHI application team also evaluate this result as the best possible calibration.

Method validation

When the NIR method was developed, there were still questions between experts and the authorities about the implementation and validation of NIR technology. For the method validation the correct identification of further 12 batches was checked. Non Topiramate tablets were also checked and an intermediate precision test performed. The robustness was demonstrated by trials with tablets that were not correctly positioned in the sample holder, with stressed tablets and placebo tablets.

Method Transfer to the new Instrument Generation NIRFlex N-500

The NIR-method "Identity of Topiramate tablets" was successfully introduced in the routine analysis at the Quality Control of Cilag AG. After the introduction, the next generation of NIR instruments and software was launched by BUCHI. The immense technical advances and the easy to use software NIRCal Toolbox signified that a change to the NIRFlex N-500 with solid transmittance module, NIRCal 5 and NIRWare Management Console software was imperative. The calibration spectra were originally measured by NIRTab instrument with a resolution of 25 cm⁻¹. For the NIRFlex N-500 with 8 cm⁻¹ resolution it was therefore necessary to convert the spectra. BUCHI undertook the data migration, imported and tested the application on the NIRFlex N-500. A validation showed that the method transfer to the latest generation of instruments and software was without any change successful and valid.

Method Robustness – Experience from the Routine

The performance of the NIR technology now became unexpectedly apparent. Between the NIR spectra of those tablets that were formerly used to develop the calibration and the batches currently manufactured, there was a slight shift. Possible reasons were different suppliers for the excipients, modified particle




Fig. 5: Report printout

sizes or tablet hardness's or different process parameters. These chemo-physical parameters that affect the NIR spectrum were covered in the original calibration. To check the effect of the shift, spectra of approximately 70 batches in all four strengths were measured (Fig. 5). All samples were positively correctly identified: the Topiramate NIR identification demonstrates to be very robust.

Calibration update

The term "ongoing method evaluation (USP)" is understood to check the performance of an existing calibration,

and is considered a possible mean of expanding the calibration with the additional incorporation of further calibration spectra. The spectra measured during routine analysis are simply allocated to the C-set (calibration spectra set) or V-set (validation set) and the calibration is recalculated. If the elimination of "older" spectra can be justified, often a higher Q-value is obtained with better performance. The updated calibration is provided to the routine user as a new application version. The application "lives"; it is also subject to a lifecycle concept in the NIRCAl software. In this way the NIR calibration is updated as per USP <1119> and Ph.EUR. 2.2.40.

Summary

The identification of the active Topiramate tablet and the discrimination between four strengths 25, 50, 100 and 200 mg is performed using a NIR transmission measurement. The time required for the IR reference method is 15 minutes, while the NIR identification requires less than 1 minute. The required calibration transfer to the latest NIRFlex N-500 model was straightforward without problems. The routine analyses have also been error free up to now and tablet batches over a period of 10 years were correctly identified. The basis for this robust identification is a calibration that was fully optimized during the development of the method. Thanks to the modular design, the transmission module for the NIRFlex N-500 basic instrument can be quickly replaced with modules for powder or liquids to use it for the identification of other, e.g. raw materials. NIR technology has come of age. The undoubtedly substantial effort required to develop the method calibration is, however, manageable for qualitative tasks. The user friendly software, good training and appropriate time resources for the laboratory analyst will remove the stumbling blocks to a successful introduction of NIR in the routine laboratory. Fascinating high technology is demonstrated to the visitor, auditor, analyst and head of laboratory: "Place sample – Scan – Ready." Now you know how Dr. McCoy in the USS Enterprise examined his patients!

BÜCHI Labortechnik AG

CH – 9230 Flawil 1
T +41 71 394 63 63
F +41 71 394 65 65
buchi@buchi.com
www.buchi.com

BÜCHI Labortechnik GmbH

DE – 45127 Essen
Freecall 0800 414 0 414
T +49 201 747 490
F +49 201 237 082
deutschland@buchi.com
www.buechigmbh.de

BÜCHI Sarl

FR – 94656 Rungis Cedex
T +33 1 56 70 62 50
F +33 1 46 86 00 31
france@buchi.com
www.buchi.fr

BUCHI UK Ltd.

GB – Oldham OL9 9QL
T +44 161 633 1000
F +44 161 633 1007
uk@buchi.com
www.buchi.co.uk

BÜCHI Labortechnik GmbH

NL – 3342 GT
Hendrik-Ido-Ambacht
T +31 78 684 94 29
F +31 78 684 94 30
netherlands@buchi.com
www.buchi.nl

BÜCHI Italia s.r.l.

IT – 20090 Assago (MI)
T +39 02 824 50 11
F +39 02 57 51 28 55
italia@buchi.com
www.buchi.it

BUCHI Hong Kong Ltd.

HK – Central
T +852 2389 2772
F +852 2389 2774
china@buchi.com
www.buchi.com.cn

BUCHI Shanghai

CN – 200052 Shanghai
T +86 21 6280 3366
F +86 21 5230 8821
china@buchi.com
www.buchi.com.cn

BUCHI (Thailand) Ltd.

TH – Bangkok 10600
T +66 2 862 08 51
F +66 2 862 08 54
bacc@buchi.com
www.buchi.com

Nihon BUCHI K.K.

JP – Tokyo 110-0008
T +81 3 3821 4777
F +81 3 3821 4555
nihon@buchi.com
www.nihon-buchi.jp

BUCHI India Private Ltd.

IN – Mumbai 400 055
T +91 22 667 18983 / 84 / 85
F +91 22 667 18986
india@buchi.com
www.buchi.in

PT. BUCHI Indonesia

ID – Tangerang 15321
T +62 21 537 62 16
F +62 21 537 62 17
indonesia@buchi.com
www.buchi.com

BUCHI Korea Inc

KR – Seoul 153-782
T +82 2 6718 7500
F +82 2 6718 7599
KR-Sales@buchi.com
www.buchi.com

BUCHI Canada Ltd.

CA – Thornhill, Ontario L4J 6Z2
T +1 416 277 7407
F +1 905 764 5139
canada@buchi.com
www.mybuchi.com

BUCHI Corporation

US – New Castle,
Delaware 19720
Toll Free: +1 877 692 8244
T +1 302 652 3000
F +1 302 652 8777
us-sales@buchi.com
www.mybuchi.com

We are represented by more than 100 distribution partners worldwide. Find your local representative at

www.buchi.com

Quality in your hands

