



Kjeldahl Weighing Boats

Reliable nitrogen/protein analysis

Weighing Boats for powder, solids or viscous fluids for reliable nitrogen and protein analysis according to Kjeldahl.

Your most important benefits:

- Convenient handling and easy transfer into sample tube
- Reliable results thanks to high quality, nitrogen-free paper



Item	Quantity	Order no.
Kjeldahl Weighing Boats, nitrogen-free paper, 58x10x10 mm	100 pieces	11060522

Contact your BUCHI partner or visit us on:
www.buchi.com/kjeldahl/consumables



www.buchi.com

Quality in your hands





Kjeldahl Tablets

Reliable nitrogen/protein analysis

BUCHI offers not only high quality instrumentation but also Kjeldahl Tablets that are tested for the most common reference applications.

Your most important benefits:

- Environmentally friendly: free of toxic mercury and selenium
- Different tablets for individual needs



Item*	Composition	Weight	Order no.
Titanium	3.5 g K_2SO_4 / 0.105 g $CuSO_4 \cdot 5 H_2O$ / 0.105 g TiO_2	3.71 g	11057980
Missouri	4.98 g K_2SO_4 / 0.02 g $CuSO_4 \cdot 5 H_2O$	5 g	11057982
ECO	3.998 g K_2SO_4 / 0.002 g $CuSO_4 \cdot 5 H_2O$	4 g	11057983
Antifoam	0.97 g Na_2SO_4 / 0.03 g Silicon Antifoam	1 g	11057984
Titanium Micro	1.5 g K_2SO_4 / 0.045 g $CuSO_4 \cdot 5 H_2O$ / 0.045 g TiO_2	1.59 g	11057981
Copper Micro	1.5 g K_2SO_4 / 0.15 $CuSO_4 \cdot 5 H_2O$	1.65 g	11057985

* each box contains 1000 tablets