

DIN/ISO/EEC/EC Methods



Nitrogen Determination according to Kjeldahl

Nitrogen Content of	DIN/ISO Method
Nitrogen content of Agricultural food products	ISO 1871
Ammonium sulphate for industrial use	ISO 3332
Milk	DIN 8968-1
Determination of non-protein-nitrogen content (NPN)	ISO 8968-4
Meat and Meat products	ISO 937
Starches & derived products	ISO 3188
Water quality, Kjeldahl Nitrogen	ISO 5663/DIN 38409H11
Water quality, ammonium	ISO 5664
Water quality, ammonium, Potentiometric method	ISO 6778
Determination of ammonia-nitrogen in water, waste water and sludge	DIN 38406-E5-2
Animal feeding stuffs - Determination of nitrogen content and calculation of crude protein content	ISO 5983-2
Rubbers Rubber, raw natural and rubber latex, natural – Determination nitrogen content	ISO 1656 (1996 E)

DIN/ISO/EEC/EC Methods



Non-Kjeldahl applications

Parameter	Method number
COD in water	DIN 38409
TVB-N in fish	EC 95/149
Phenol determination in sludge and beer	ISO 6439: DIN 38409-H16-3
Alcohol determination in wine, sprits	EEC 2676/2870