



# Vacuum Pump V-710

## The powerful vacuum source

The V-710 is a powerful, chemical resistant and compact vacuum pump. The silent system is optimally suited for Industrial Rotavapor® or to act as a central lab vacuum source for multiple Rotavapor®.



### Powerful

Fast and low vacuum generation thanks to the powerful four diaphragm heads technology



### Flexible

Covers a wide range of different applications



### Robust

Long-lasting and chemical resistant due to exclusive use of high quality materials

## V-710: Your most important benefits



### Powerful

- Powerful pump with four diaphragms for multiple use in your lab
- Fast evacuation of even high volume Rotavapor® due to the high suction volume



### Flexible

- Seamless integration into the BUCHI Rotavapor® system
- Adequate alternative to central vacuum source
- Easy control and best performance in combination with a BUCHI Vacuum Controller



### Robust

- Robust and chemical resistant construction for best longevity
- Nearly maintenance-free due to exclusive user of high quality materials
- Increases life time thanks to integrated ECO<sup>2</sup> mode

## Complete your portfolio



Vacuum Controller  
V-850 / V-855  
Achieve best process  
performance



Rotavapor®  
R-220 SE  
Best choice for  
your scale up lab



Rotavapor®  
R-215  
Excellence in  
evaporation



Multivapor  
P-6 / P-12  
Efficient evaporation  
for multiple samples

