



Syncore® Analyst

Economical high quality concentration

The Syncore® Analyst offers efficient, gentle and environmentally friendly parallel concentration. With its clever design cross-contamination can be inhibited and recoveries are maximized. Profit from the wide range of complementary products to tailor your sample workup to your needs.

Efficient

Increase productivity by concentrating up to 12 samples simultaneously



High quality results

Maximize recoveries and eliminate the risk of cross contamination



Time savings

Fast Vortex concentration maximizes throughput to a predefined residual volume

Syncore® Analyst: Your most important benefits



Efficient

- Improve your sample throughput by concentrating up to 12 samples to a residual volume
- No need of nitrogen saves operational cost
- Small foot print saves expensive lab space
- Optimize your process by combining the solid phase extraction and the concentration



High quality results

- Eliminate the risk of cross contamination due to a clever design
- Integrated solid phase extraction (SPE) eliminates contamination risks
- High recovery rates due to smooth sample processing (up to 100 % solvent recovery)



Time savings

- Save time by using the fast and proven parallel Vortex evaporation technology
- Increase productivity by fast concentration of up to 12 samples in parallel
- Concentrate to a predefined residual volume allows a unattended workup

Complete your portfolio



SpeedExtractor
E-916 / E-914
Pressurised solvent
extraction



Syncore®
SPE
Parallel concentration and
solid phase extraction



Vacuum Pump
V-300
The economical and
silent vacuum source



Recirculating Chiller
F-305 / F-308 / F-314
The efficient and water
saving way of cooling

